

dPette

Simple Electronic Pipette

Features

- Motor driven digitally control pipette with multifunctions
- Easy Operations and Handling
- Intuitive menu interface settings: functions and operations
- Minimal force for pipetting operation
- High performance ensuring accuracy & repeatability
- 2 buttons for all operational settings
- Adjustable speed for aspiration and dispensing
- Li-ion battery enable longer operation time
- Convenient and Versatile
- Self-calibration applicable



dPette
For Pipetting
and Mixing

Specifications

Channels	Volume Range	Increment	Test Volume	Systematic Error		Random Error	
	μL	μL	μL	μL	%	μL	%
1	0.5-10	0.01	10	± 0.10	± 1.00	± 0.05	± 0.50
			5	± 0.10	± 2.00	± 0.10	± 2.00
			1	± 0.05	± 5.00	± 0.03	± 3.00
1	5-50	0.1	50	± 0.40	± 0.80	± 0.15	± 0.30
			25	± 0.25	± 1.00	± 0.25	± 1.00
			5	± 0.20	± 4.00	± 0.125	± 2.50
1	30-300	1	300	± 1.80	± 0.60	± 0.60	± 0.20
			100	± 1.00	± 1.00	± 0.40	± 0.40
			30	± 1.20	± 4.00	± 0.21	± 0.70
1	100-1000	5	1000	± 6.00	± 0.60	± 2.00	± 0.20
			500	± 5.00	± 1.00	± 1.00	± 0.20
			100	± 3.00	± 3.00	± 0.60	± 0.60

* DLAB specifications are used as guidelines and the user calibration should refer to the industrial standard ISO 8655