

Features

- Available as Fixed Volume and Adjustable Volume formats
- Streamlined shape design
- Comprehensive volume range from 0.1µL to 10mL
- Easy calibration and maintenance
- 8 and 12 channel pipette options
- Dispensing head rotates for effortless pipetting convenience
- Individual piston and tip cone assembly
- Spring loaded tip cones for easy cleaning and maintenance
- Compound material-made tip cone secures high sealing performance



8-channel Adjustable Volume Pipettes

Volume Range	Increment	Test Volume	Systematic Error		Random Error	
			µL	%	µL	%
0.5-10	0.1	10	±0.15	±1.50	±0.15	±1.50
		5	±0.125	±2.50	±0.125	±2.50
		1	±0.04	±4.00	±0.04	±4.00
5-50	0.5	50	±0.5	±1.00	±0.25	±0.50
		25	±0.375	±1.50	±0.25	±1.00
		5	±0.15	±3.00	±0.1	±2.00
10-100	1	100	±0.8	±0.80	±0.3	±0.30
		50	±0.5	±1.00	±0.4	±0.80
		10	±0.3	±3.00	±0.2	±2.00
50-300	5	300	±2.1	±0.70	±0.75	±0.25
		150	±1.5	±1.00	±0.75	±0.50
		50	±0.75	±1.50	±0.4	±0.80
100-1000	10	1000	±10	±1.00	±3	±0.30
		500	±10	±2.00	±2.5	±0.50
		100	±5	±5.00	±1	±1.00

12-channel Adjustable Volume Pipettes

Volume Range	Increment	Test Volume	Systematic Error		Random Error	
			µL	%	µL	%
0.5-10	0.1	10	±0.15	±1.50	±0.15	±1.50
		5	±0.125	±2.50	±0.125	±2.50
		1	±0.04	±4.00	±0.04	±4.00
5-50	0.5	50	±0.5	±1.00	±0.25	±0.50
		25	±0.375	±1.50	±0.25	±1.00
		5	±0.15	±3.00	±0.1	±2.00
10-100	1	100	±0.8	±0.80	±0.3	±0.30
		50	±0.5	±1.00	±0.4	±0.80
		10	±0.3	±3.00	±0.2	±2.00
50-300	5	300	±2.1	±0.70	±0.75	±0.25
		150	±1.5	±1.00	±0.75	±0.50
		50	±0.75	±1.50	±0.4	±0.80