

# SK-0180-Pro

## LCD Digital Orbital Shaker

### Features

- Independent LCD display for speed and time
- Compact footprint
- Provides smooth and uniform shaking
- Two options: 7.5 kg (O330) with platform and 2.5kg (O180) with platform speed range 100-500 rpm and 100-800 rpm
- Continuous or timer operation mode
- Maintenance- free brushless DC motor
- Over speed detection and protection
- Wide range of optional platform for flasks, and other popular vessels
- Remote function provides PC control and data transmission
- Free software for PC control & data logging



**SK-0180-Pro**

### Applications











Suitable for mixing in culture dishes, flasks and beakers. Ideal for aeration and extraction applications. It can be used in environmental chambers and CO2 incubators.












Orbital Motion



Specifications	SK-O180-Pro
Orbital diameter	4mm
Max. load capacity (with platform)	2.5kg
Motor type	Brushless DC motor
Motor input power	28W
Motor output power	15W
Speed range	100-800rpm
Speed display	LCD
Timer display	LCD
Time settings range	1min-99h59min
Voltage,Frequency	100-240V,50/60Hz
Power	30W
Weight	8.1kg
Permissible ambient temperature and humidity	5-40°C,80%RH
Dimension[W×D×H]	300×340×100mm
Protection class	IP20
Data connector	RS232

Cat. No.	Description	SK-O330-Pro	SK-L330-Pro	SK-R330-Pro	SK-D3309-Pro	SK-O330-M	SK-O180-Pro	SK-L180-Pro	SK-O180-S Multifunctional models
 18900036	Bar for SK330.1 universal attachment	■	■	■	■				
 18900037	Bar for SK180.1 universal attachment						■	■	■
 18900280	Connect sticks for double dishes shaker	■	■	■	■	■	■	■	■
 18900028	SK330.2, Fixing clip attachment Size: 32 x 32cm	■	■	■	■				
 18900026	SK180.2, Fixing clip attachment Size: 25 x 24cm						■	■	■
 18900029	Fixing clip for flask volume 25mL, used with SK330.2[9pcs] SK180.2[28pcs]	■	■	■	■		■	■	■
 18900030	Fixing clip for flask volume 50mL, used with SK330.2[9pcs] SK180.2[16pcs]	■	■	■	■		■	■	■
 18900031	Fixing clip for flask volume 100mL, used with SK330.2[9pcs] SK180.2 [16pcs]	■	■	■	■		■	■	■
 18900032	Fixing clip for flask volume 200/250mL, used with SK330.2[4pcs] SK180.2[9pcs]	■	■	■	■		■	■	■
 18900033	Fixing clip for flask volume 500mL, used with SK330.2[2pcs] SK180.2[9pcs]	■	■	■	■		■	■	■

Cat. No.	Description	SK-O330-Pro	SK-L330-Pro	SK-R330-Pro	SK-D3309-Pro	SK-O330-M	SK-O180-Pro	SK-L180-Pro	SK-O180-S Multifunctional models
 18900155	SK330.5, Dish platform with non-slip mat and fixed string (8pcs) Dimension[WxD]: 30x28.5cm	■	■	■	■				
 18202835	Silicone mattress for 16× 5~15mL tubes, use with SK330.5	■	■	■	■				
 18202846	Silicone mattress for 6× 50mL tubes, use with SK330.5	■	■	■	■				
 18900185	SK180.4, Dish attachment with non-slip mat Dimension[WxD]: 24x21.5cm						■	■	■
 18900869	SK330.4 Dish platform with non-slip mat and fixed string(8 pcs), 34X34cm					■			
 18900859	SK330.6 Transparent shield, serologic test for syphilis, 41X39.5X10cm	■	■			■			
 18900860	SK180.5 Transparent Cover for Syphilis Serum Test, 32X32X8.3cm						■	■	
 18202852	Silicone mattress for 10× 5~10mL tubes						■	■	■
 18202853	Silicone mattress for 4× 50mL tubes						■	■	■
 18900027	SK330.1, Universal attachment with 4 bars Size: 32x32cm Capacity for flask: 25mLx21pcs, 50mLx15pcs, 100mLx12pcs, 250mLx6pcs, 500x4pcs	■	■	■	■	■			
 18900040	SK330.3 Lengthways roller attachment, Mainly used for separation funnel oscillation extraction applications 500mL*3pcs 250mL*4pcs 100/125mL*4pcs 50mL*5pcs	■	■	■	■	■			
 18900025	SK180.1, Universal attachment with 3 bars Size: 25x24.5cm Capacity for flask: 25mLx10pcs, 50mLx8pcs, 100mLx6pcs, 250mLx4pcs, 500x2pcs						■	■	■
 18900038	SK180.3 Lengthways roller attachment, Mainly used for separation funnel oscillation extraction applications 500mL*2pcs 250mL*3pcs 100/125mL*3pcs 50mL*4pcs						■	■	■